

Standard Operating Procedure (SOP) 005V7.0

Processing and Storage of Breast Tissue

SPREC TIS-BPS-N-B-SNP-A-C [1]

SPREC TIS-BPS-N-B-PXT-G-P [1]

Effective Date: March 2, 2020

Author: Pamela Rockey, M.S., Biospecimen Manager

Approved by:



Anna Maria Storniolo, M.D., KTB Executive Director

Materials:

Cryostorage tubes: Corning 2.0ml Cryogenic Vials. (Fisher cat. #0337421)

Bleach: Any commercially available

Polypropylene Specimen Containers: (Fisher cat. #14 955 103)

Liquid Nitrogen: (ICTSI-SSF R3 C156)

Forceps: (Fisher cat. #08-887)

Embedding (Tissue) Cassettes: (Fisher cat. #15-182-708B)

PAXgene Tissue Fix Containers: (Qiagen #765312)

PAXgene Tissue Stabilizer: (Qiagen #765512)

Ethanol: Molecular Biology Grade Absolute

Stiff Bristled Brushes: (Fisher #03-635)

CryoPod™ Carrier (Brooks Life Sciences)

Balance and Weigh Papers

Emesis Basin or a receptacle to hold diluted bleach and forceps

Appropriate PPE: Gloves, Lab Coat, Eye protection

Normal breast tissue is procured following SOP 001V8.0 (Acquisition of Normal Breast Tissue and Blood at a Tissue Collection Event) or SOP 012V6.0 (Acquisition of Normal Breast Tissue and Blood from an Operating Room) from donors meeting all eligibility criteria as listed here: <https://komentissuebank.iu.edu/donate-tissue/register-to-donate.php>

Tissue Processing:

FIXATION: Within 10 minutes (or less) of procurement tissue cores are received from the biopsy room on a telfa pad and separated into individual cores using a clean forceps. One breast tissue core is placed into an embedding cassette and the cassette is closed, weighed and weight is recorded. (Tissue cassettes are labeled in #2 pencil lead or by bar code label prior to tissue acquisition.) As described in the PAXgene Tissue FIX Container (50 ml) Product Circular [2], the cassette is placed into a 50 ml PAXgene FIX tissue

container where the tissue is fixed. The sample ID and time is written on the top label and on the appropriate spreadsheet. Four (4) cassettes are added to each 50 ml FIX container.

STABILIZATION: Two hours after the 4th or last cassette has been placed in the container, the FIX solution is poured into a chemical waste container under a chemical fume hood. Care is taken not to let the cassettes fall out during pouring. Approximately 50ml of PAXgene Stabilizer (prepared from concentrate) is poured into the container ensuring tissues are covered and the lid is secured. The time of Stabilizer exchange is recorded on the side of the container and on the appropriate spreadsheet. The tissue container is sealed in its original plastic bag and allowed to remain at room temperature. At the end of the collection event, tissue containers are transported to the KTB lab and allowed to remain at room temperature for up to 72 hours. If the cups have not been sent to processing within 72 hours, they can be stored at 2-8°C for up to 4 weeks. Tissue containers are shipped to a processing facility where the PAXgene processing protocol for paraffin embedding is strictly adhered to. A 5µm section of the resulting paraffin block is obtained for hematoxylin and eosin staining. A digital image of the H&E is entered into the database.

Universal Precautions are MANDATORY. Eye protection is mandatory every time liquid nitrogen is handled to protect against injury due to splashing. Standard laboratory personal protective equipment (e.g. closed toe shoes, full cover of legs and feet, and goggles) will be worn when handling coolants.

FRESH CORES: If requested by investigator, **within 10 minutes (or less) of procurement** one or two breast tissue cores are placed in investigator supplied media in labeled 50 ml conical tubes and transported immediately at RT or on wet ice to the on-site processing lab.

FRESH FROZEN: Flash freezing in liquid nitrogen provides excellent specimen integrity and a wide array of options for tissue analysis [3]. **Within 10 minutes (or less) of procurement**, enough liquid nitrogen (LN₂) to cover the specimen is transferred to a new polypropylene container. Using a clean forceps, the remaining cores are immediately placed into LN₂ and snap frozen for 30-60 seconds. The frozen cores are placed in labeled, chilled cryovials. (All cryovials pre-labeled with bar-coded labels prior to tissue acquisition). The cryovials are scanned with a barcode reader, logged into cryoboxes and held in a charged CryoPod™ Carrier or charged LN₂ dry shipper for transport to the storage facility.

Used forceps are placed in a bleach solution (1/10 bleach/water) in an emesis basin. When time allows or a clean forceps is needed, the forceps in the bleach solution are gently scrubbed with a brush, rinsed in water and allowed to air dry.

Temperature of specimens following acquisition and Snap Freezing but prior to storage: Cryovials containing the specimens are to be held at ≤ -150°C in an LN₂ CryoPod™ Carrier or LN₂ dry shipper until transfer to final storage location.

Storage of Tissue: Once all the samples are procured, the cryovials are transported in CryoPod™ Carriers or LN₂ dry shipper to the SSF LN₂ tank. Tissue samples are stored in liquid nitrogen vapor, (-166.2°C to -195.1°C from top to bottom of tower).

Temperature for collection and processing: All tissue procurement and processing events are done at room temperature.

Standardization: All variables including the time between excision and snap-freezing and time stored in LN2 vapor phase prior to shipment and/or utilization will be entered into the database.

Oversight: All adverse and unexpected events will be recorded in the database and will be addressed by the Internal Advisory Committee. This includes all phases of the process: donation, processing, storage and retrieval, and utilization.

References:

1. Sabine Lehmann et.al. International Society for Biological and Environmental Repositories (ISBER) Working Group on Biospecimen Science. Standard Preanalytical Coding for Biospecimens: Review and Implementation of the Sample PREanalytical Code (SPREC). *Biopreservation and Biobanking* Vol. 10 No.4, 2012
2. PAXgene® Tissue FIX Container (50ml) Product Circular. Feb 2013. PreAnalytiX a Qiagen/BD Company. www.preanalytix.com
3. https://biospecimens.cancer.gov/global/pdfs/NCI_BEBP_Snap-freezing_of_Post-surgical_Tissue_Biospecimens.pdf

Bibliography

- Karlsson JO, Toner M. *Long-term storage of tissues by cryopreservation: critical issues*. *Biomaterials*. 1996; **17**(3):243-56
- Morente MM, Mager R, Alonso S, Pezzella F, Spatz A, Knox K, Kerr D, Dinjens WN, Oosterhuis JW, Lam KH, Oomen MH, van Damme B, van de Vijver M, van Boven H, Kerjaschki D, Pammer J, Lopez-Guerrero JA, Lombart Bosch A, Carbone A, Gloghini A, Teodorovic I, Isabelle M, Passiukov A, Lejeune S, Therasse P, van Veen EB, Ratcliffe C, Riegman PH. *TuBaFrost 2: Standardising tissue collection and quality control procedures for a European virtual frozen tissue bank network*. *Eur J Cancer*. 2006; **42**(16):2684-91.
- Qualman SJ, France M, Grizzle WE, LiVolsi VA, Moskaluk CA, Ramirez NC, Washington MK. *Establishing a tumour bank: banking, informatics and ethics*. *Br J Cancer*. 2004; **90**(6):1115-9
- Villalba R, Eisman M, Fornes G, Mesa A, De Torres P, Ariza MJ Gutierrez Guzman A, Gomez Villagran JL. *Implementation of a Quality Plan (ISO 9002) In a Regional Tissue Bank*. *Cell Tissue Bank*. 2001; **2**(1):45-49.
- Leyland-Jones, Brian R. et al., *Recommendations for Collection and Handling of Specimens from Group Breast Cancer Clinical trials, from Onsite Collection through Shipping to the Central Bank*. *Journal of Clinical Oncology* 2008 26:34, 5638-5644
- Campbell LD, Astrin JJ, DeSouza Y, Giri, J, Patel AA, Rawley-Payne M, Rush A and Sieffert N. *The 2018 Revision of the ISBER Best Practices: Summary of Changes and the Editorial Team's Development Process*. *Biopreservation and Biobanking* 16(1): 3-6.
<https://doi.org/10.1089/bio.2018.0001>

Electronic Resources

- ISBER Best Practices: Recommendations for Repositories
<https://www.isber.org/page/BPR>
- NCI Best Practices for Biospecimen Resources:
<https://biospecimens.cancer.gov/bestpractices/2016-NCIBestPractices.pdf>
- Snap-Freezing of Post-Surgical tissue Biospecimens:
https://biospecimens.cancer.gov/global/pdfs/NCI_BEBP_Snap-freezing_of_Post-surgical_Tissue_Biospecimens.pdf